

Labycare-Tube Sealer Gun Type (SKSS388A)



Portable gun-type tube sealer for fast, reliable sealing of blood bag tubing. Features adjustable sealing time, rechargeable battery, and safety cut-off for efficient, leak-proof operation in blood banks and laboratories.

Description :

Product Cat # SKSS388A

Specifications :

Parameter

Sealing Method
Seal Width
Sealing Time
Compatible Tubing
Power Supply
Battery Type
Seals per Charge

Specification

Thermal impulse sealing
5–7 mm
Adjustable (0.5 – 2.0 seconds)
PVC and medical-grade tubing
100–240 V AC, 50/60 Hz (with adapter)
Rechargeable Li-ion battery
250–400 seals (approx.)

Charging Time	2–3 hours
Indicators	LED for Power, Ready, and Seal Status
Safety Features	Overheat protection, automatic cut-off
Dimensions	Approx. 180 × 60 × 40 mm
Weight	Approx. 350 g
Application	Blood banks, hospitals, laboratories

labycare

www.labycare.com | Email: labycare@gmail.com