

Labycare-Platelet Incubator Agitator. (SKSS389)



Platelet Incubator Agitator with precise temperature control and gentle agitation, ensuring optimal platelet storage and extended viability for transfusions.

Description :

The Platelet Incubator Agitator ensures optimal storage of platelet concentrates with controlled temperature and gentle agitation. Designed for blood banks and laboratories, it maintains platelet viability, enhances shelf life, and ensures uniform conditions, supporting safe and effective transfusions.

Specifications :

Feature	Specification
Temperature Range	20°C – 24°C ($\pm 0.5^\circ\text{C}$)
Agitation Speed	60 – 80 cycles/min
Capacity	4–6 platelet bags per shelf (varies by model)
Control	Digital LED display with alarms
Power Supply	220V \pm 10%, 50Hz
Safety Features	Over-temperature alarm, power failure alarm
Construction	Corrosion-resistant stainless steel

labycare

www.labycare.com | Email: labycare@gmail.com