

Labycare-Plasma With Solar System Refrigerator. (SKSS383C)



Based on the unique Solar chill ice bank technology, the labycare.

Description :

Based on the unique Solar chill ice bank technology, the Vest frost MKS 044D functions safely without use of energy storing batteries or contaminating fuel. Solar panels supply the necessary power and the ice bank provides a reliable temperature in the vaccine compartment- even when the sun is absent for days. The temperature is controlled and maintained via self-regulating convection between ice bank and vaccine compartment

<u>Specifications :</u> Temperature Range	+ 2 to +8 Degree
Model Name/Number	SKSS383C
Usage/Application	Laboratory
Size	L28" X W28" X H33"

labycare

www.labycare.com | Email: labycare@gmail.com