

Labycare-Blood Donor Bed Couch. (SKSS820)



We design **Blood Donor Couch** with perfection to ensure comfort of the blood donor. Our Blood Donor Couch ensures safe transfusion of blood ensuring desired laying posture of the blood donor. Our Blood Donor Couch has wheels attached to it, which ensure easy movement from one place to another. We make available uniquely designed Blood Donor Couch at market leading prices.

Description :

- Facility for blood collection from both sides
- Provides a comfortable position for the donor
- 4 Heavy Duty Lockable castors for easy mobility
- Based on haemodynamic principles. (Donor's head can be lowered immediately and legs can be lifted above heart level so that blood can flow back to The brain and other vital organs, in case of vasovagal attack)
- Variable positioning for either arm with a comfortable wide arm-rests
- Tilt adjustment can be done by using remote control
- Comfortable chair type with soft padding for cushioning and rexin cover
- Two Geared Actuator
- Ensuring safety and comfort to the donor

Specifications :

~~Copyright © 2019 Labycare. All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording, or by any information storage and retrieval system, without prior written permission from Labycare.~~

labycare

www.labycare.com | Email: labycare@gmail.com