

Labycare-Hospital Beds (SKSS801M)



Hospital Beds manufacturer by Labycare in India. It offers the greatest convenience for the patient and caregiver. The easy-to-use screw liner mechanism helps in positioning of the upper body and knees. ABS headboard and foot board are more durable, impact and scratch resistant, washable for easy cleaning. 5" Lockable casters helps in easy mobility and reduce slip. Safety side rails helps prevent patients from falling and gives them a secure feeling. 4" waterproof mattress will be against water spills.

Description :

Our rich industry experience enables us to manufacture and export an extensive range of **Labycare Hospital Bed Simple**, which are designed as per the international quality standards. Moreover, these beds are offered after stringently checking and analyzing them on well defined parameters of the quality. Our range is widely accredited for their smooth finishing, durability, comforts, adjustable nature and reliability.

Specifications :

Model NO.

SKSS801M

Condition

New

Name

Stainless Steel Two Cranks Manual Hospital Bed

Service Manual

Available

User's Manual

Available

Warranty

2 Years

Samples

Available

Certifications

-

Factory Location

Ambala cantt

Service

Lifetime

Spare Parts

1% Free

Factory Type

Over 13 Years Manufacturer

Trademark

Labycare

Transport Package

Carton/Wooden Box/Others

Specification

L2150*W960*H500 mm

Origin

India

labycare

www.labycare.com | Email: labycare@gmail.com