

## Labycare-Electrical ICU Bed 4 Functions (SKSS801)



Electrical ICU bed with four functions provides motorized adjustment of backrest, knee rest, height, and Trendelenburg positions for safe, comfortable, and efficient critical care management.

### **Description :**

An electrical ICU bed with four functions is designed for critical care environments, providing smooth motorized adjustment of backrest, knee rest, height, and Trendelenburg positions. It enhances patient comfort, supports efficient nursing care, and ensures safety and accessibility during intensive monitoring and treatment.

### **Specifications :**

Functions: Backrest, knee rest, height adjustment, Trendelenburg / Reverse Trendelenburg

Operation: Electric motor with hand control / touch panel

Frame: Heavy-duty mild steel structure

Bed surface: Perforated CRCA / ABS platform

Side rails: Collapsible ABS safety rails

Finish & mobility: Epoxy powder coated with lockable castor wheels

---

labycare

www.labycare.com | Email: labycare@gmail.com