

Labycare-Electrical ICU Bed 3 Function (SKSS801N)



Electrical ICU Bed 3 Function Manufacturer by LABYCARE

Description :

Electrical ICU Bed with 3 motorized functions is designed for critical care units, offering smooth height, backrest, and leg section adjustment. Built with a sturdy frame, easy mobility, and patient-friendly design, it ensures comfort, safety, and efficient caregiver operation in intensive care environments.

Specifications :

- Three electric functions: backrest, knee rest, and height adjustment
- High-strength mild steel frame with epoxy powder coating
- Four-section mattress platform with radiolucent backrest
- Lockable castor wheels with smooth mobility
- Detachable ABS head and foot panels
- Manual CPR release for emergency use

labycare

www.labycare.com | Email: labycare@gmail.com