

Labycare-Full Fowler Bed With ABS Panel (SKSS809C)



Hospital Semi Fowler Beds with ABS Panel manufacturer by **Labycare** in India. We are leading manufacturer of Full fowler bed with abs panel. We can provide you best quality product at very reasonable rate.

Description :

Specifications :

- Mild steel tubular Head and Foot Bows of unequal height with vertical tube support
- Framework made up of M.S. tubes
- 4 section top made up of oval perforated CRCA M.S. sheet
- Bed mounted on 4" Castors, 2 with brakes
- Four location for I.V. Rod
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2000L x 910W x 600H mm
- Positions obtained by Smooth Crank System with foldable handle:

labycare

www.labycare.com | Email: labycare@gmail.com