

Labycare-Semi Fowler Bed With ABS Panel (SKSS809D)



Labycare Semi Fowler Beds with ABS Panel manufacturer by **Labycare** in India. We are the one of best and foremost manufacturer of Semi Fowler bed in India. We do customize the product as per client request or need. We do use the best quality of materials in semi Fowler bed for benefit of client and patient. Semi Fowler bed is an important feature to be used in hospital stay of patient. So we do our best to help you in fast recovery of your patient.

Description :

SK Scientific & Surgicals is considered as a renowned manufacturer and supplier of Hospital Bed.

Specifications :

- ABS moulded Head and Foot Panel with locking and corner buffer
- Framework made up of M.S. tubes
- 2 section top made up of oval perforated CRCA M.S. sheet
- Legs fitted with PVC stumps
- Four location for I.V. Rod
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2000L x 910W x 600H mm
- Positions obtained by Smooth Crank System with foldable handle:

1. Backrest elevation

Optional

- 2 Section Mattress
- S.S. Telescopic I.V. Rod
- S.S. Collapsible Railing
- Advance linear actuator system (01 Pc Motor, 01 Pc control box, 01 Handset)
- M.S. or S.S. Head & Foot Panels

labycare

www.labycare.com | Email: labycare@gmail.com