

Labycare-Hospital Semi Fowler Bed (SKSS809Q)



Hospital Semi Fowler Beds manufacturer by Labycare in India.

Description :

Specifications :

Bed type: Semi Fowler hospital bed

Operation: Manual, crank-operated

Adjustable sections: Backrest and knee rest

Frame: Mild steel tubular construction

Bed surface: Perforated CRCA sheet

Finish: Epoxy powder coated

Side rails: Optional / collapsible

Mattress size: Approx. 1980 × 900 mm

Load capacity: Approx. 150–180 kg

Application: General wards and nursing care units

