

Labycare-Hospital Plain Bed (SKSS809E)



Hospital Plain Bed manufacturer by **Labycare** in India. Owing to the efforts of our dedicated team of professionals, we have been constantly engaged in offering the best quality **Hospital Plain Bed**.

Description :

We offer an extensive range of **Labycare Hospital Plain Bed**. Our beds are strong due to the precise steel tubular frame work. These beds are mostly used in hospitals, nursing homes & health centers and are acknowledged for excellent durability. It has tubular Head and Foot Bows of unequal height. It comes with Pretreated and Epoxy powder coated finish and is available at very competitive prices.

Specifications :

- Mild steel tubular Head and Foot Bows of unequal height with vertical tube support
- Framework made up of M.S. tubes
- Top made up of oval perforated CRCA M.S. sheet
- Legs fitted with PVC stumps
- Four location for I.V. Rod
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 1900L x 910W x 600H mm

- Plain Mattress
- S.S. Telescopic I.V. Rod
- S.S. Collapsible Railing

labycare

www.labycare.com | Email: labycare@gmail.com