

Labycare-Hospital Labour Table (SKSS818D)



Hospital Labour Table manufacturer by **Labycare** in India.

Description :

We are a company which is well known in the market for offering a wide range of **Labycare Labour Table** which is used in hospitals, medical centers and other allied places to provide comfort to the patient having labor pain. It is made available with easily adjustable backrest and other required assemblies. This range is designed for use in maternity wards and labour rooms. Offered **Labour Table** is appreciated for attractive design, Adjustable Trendelenburg positions, strong construction and supreme finish.

Specifications :

- Overall approx. size: 72" x 27" x 30"
- Frame work made up of S.S. tube
- Legs fitted with PVC stumps
- 3 section top made up of Bakelite Sheet
- Back rest section adjustable ratchet mechanism from back end
- Leg end section can slide under the main section
- The entire table top gives trendelenburg positions by gas spring from back end
- Drop down S.S. safety side railing are provided on 3 sides of the table
- Height adjustable a pair of foam padded knee crutches
- Patient hand grips made up of S.S. tubes
- Supplied with waste collecting S.S. Bowl on S.S. Bracket and one S.S. I.V. Rod
- Also available in M.S. Frame

labycare

www.labycare.com | Email: labycare@gmail.com