

Labycare-Delivery Telescopic Table (SKSS818C)



Delivery Telescopic Table manufacturer by **Labycare** in India. To accomplish the specific demands of our clients, offered table is available in varied specifications. This **Telescopic Labour Table** can be availed from us at economical prices.

Description :

Labycare Telescopic Labour Table is extensively demanded in hospitals, gynecology clinics and other allied places to lay patients for delivery purpose. Our provided table is manufactured by utilizing high grade raw materials with the aid of latest technology in accordance with global industry standards. To accomplish the specific demands of our clients, offered table is available in varied specifications. This **Telescopic Labour Table** can be availed from us at economical prices.

Specifications :

- Approx. dimensions: 72"L x 27"W x 24"-32"H
- Tubular frame work mounted on swiveling castors
- Bakelite sheet top in three sections with 3 section mattress
- Trendelenburg / Reverse Trendelenburg positions, adjustable by gas spring from back end
- Hydraulically height adjustable by foot paddle from side
- Backrest section adjustable crank mechanism from back end
- Leg end section can slide under the main section
- Height adjustable by a pair of knee crutches
- Pre-treated & Epoxy Powder Coated
- Also available in S.S. Frame & S.S. Top