

Labycare-Gynecological Examination Couch (SKSS822)



Gynecological Examination Table manufacturer by **Labycare** in India.

Description :

The Gynecological Examination Table is designed for comfort, functionality, and ease of use during gynecological and obstetric examinations. Made with a robust steel frame, cushioned upholstery, and adjustable sections, it ensures both patient comfort and practitioner convenience

Specifications :

- **Material:** Mild steel frame with epoxy powder coating
- **Dimensions:** 1800 x 600 x 800 mm (L x W x H)
- **Upholstery:** High-density foam with washable, antimicrobial cover
- **Adjustable Features:** Backrest and leg support with smooth mechanical adjustment
- **Leg Supports:** Padded knee crutches with height and angle adjustment
- **Trendelenburg Position:** Available with manual adjustment
- **Accessories:** Optional stainless steel tray and lithotomy rod holders
- **Weight Capacity:** 150 kg
- **Customization:** Available for dimensions and additional features.

labycare

www.labycare.com | Email: labycare@gmail.com