

## Labycare-Hospital Dressing Trolley (SKSS830)



**Hospital Dressing Trolley** manufacturer by **Labycare** in India and having a definite quality management system, we are increasingly engaged in providing an inclusive **Hospital Dressing Trolley**.

### Description :

**Labycare Hospital Dressing Trolley** is mainly used in hospitals / medical centers to carry medicines and surgical instruments. These trolleys are generally square metal trolley with rollers on each leg to assist in the quick movement of the trolley around the place. It also used for **dressing surgical procedures in wards**. Stainless Steel Trolley with Rails. Dressing Trolleys consist of stainless steel or epoxy coated mild steel frames, and usually comes with two stainless steel shelves, and a rail. In line with the ever growing desires of our patrons, we have come up with a comprehensive consignment of Surgical Dressing Trolley. Obtainable with us in a consignment of sizes and dimensional stipulations, these presented dressing trolleys are broadly recommended and acknowledged. Together with this, our patrons can avail these from us in a consignment provisions at affordable costs.

### Specifications :

- Approx. size : 30" x 18" x 32"
- Frame work made up of S.S. tubes
- Two S.S. shelves up with S.S. railings on all four sides
- The trolley fitted S.S. bucket and S.S. bowl
- The trolley mounted on 100 mm dia castors
- S.S. frame, S.S. shelves with S.S. basin & bucket
- Also available in Epoxy Powder Coated Frame & S.S Top

---

labycare

www.labycare.com | Email: labycare@gmail.com