

Labycare-Hospital Over Bed Table (SKSS837)



Hospital Over Bed Table is manufacturer by **Labycare** in India. It is easy to use for any Hospital Bed. Height adjustable table top, requires minor assembly.

Description :

- Chrome-plated, welded tubular steel column and base provide durability and stability.
- 1.5" hooded, ball-bearing casters ensure easy maneuverability across all floor surfaces.
- Table top, constructed of wood-grained oak Melamine laminate on particle-board core, adjusts 13" in infinite positions.
- Raised T-molded plastic edge.
- Table width: 30"
- Table Depth: 15"
- Height adjustment: 30" - 43"
- Chrome-plated "H" base and column for aesthetics and durability.
- 2 locking casters
- Supports up to 30 lb of dispersed weight.

Specifications :

Material Of Top	Mild Steel
Length Of table top	900 millimeter
Width Of table top	600 millimeter
Thickness Of table top	1 millimeter
Design of Over bed Table	Ability to position bed table top over patient from side of bed
Provision for changing slope of table top from its horizontal position	Yes
Material Of Frame	Mild Steel
Cross Section Of Frame Material	Square

Safe Working Load (maximum)	15 Kg or more.
Dia/size of pipe used in frame	16 millimeter
Thickness Of Pipe Of Frame	1.2 millimeter

labycare

www.labycare.com | Email: labycare@gmail.com