

Labycare-Modular Cold Room (SKSS374)



Modular cold rooms manufacturer by **Labycare** in India. They can be partitioned or coupled and are available in varying thicknesses and heights.

Description :

Labycare Modular Cold Rooms are Insulated chambers built at site to desired dimensions and cooling specifications to meet our customer's requirement. Cold rooms provide precise temperature control for commercial facilities that need powerful, reliable and consistent refrigeration or freezing – and for unstable perishable items like food and chemicals, this means long-term temperature control is needed.

Specifications :

- **Outside Dimension:** 8'-0"Front x 8'-0"Side
- **Finish:** 26 gauge stucco galvanized steel
- **Wall Panels:** Metal-clad, framed panels with cam lock fasteners, cam lock walls, ceilings and floors.
- **Insulation:** 4 inch foamed in place urethane with a R-28 value.
- **Floor not Included Door (s)** are 36" wide x 78" high, self closing with magnetic gasket seal and with spring loaded cam lift chrome hinges also includes.
- **Walk- in Door** a locking handle. 1 door viewport not included, interior kickplate not included, and Exterior kickplate not included.
- **Front** one 2-1/2 inch dial thermometer. Two vapour proof LED light fixture. Two pilot light switches.

labycare

www.labycare.com | Email: labycare@gmail.com