

## Labycare-General Operating Table (SKSS859)



- **General Operating Table** manufacturer by **Labycare** in India. Table's design gives a surgical team the flexibility, control and convenience they need for today's patients and procedures. Provides outstanding access for C-arm, permitting high quality images crucial to cardiothoracic, neurosurgical and orthopedic procedures. Facilitates quick and safe patient positioning with an easy-to-use, easy-to-read hand control contours to virtually all patient postures with five section tabletop.

### **Description :**

**Labycare General Operating Table** is high quality operating table with built in multiple functions making it suitable for general surgery Functional design & easy to use hydraulic adjustment Wide accessory range for specialized procedures 3 section table top frame of Stainless steel Antistatic PU Integral mattress Height around 730 – 1040 mm. Trendelenburg 25 — +25 degree Side tilt +/- 15 deg with gas spring Back section -15 – 70deg with gas spring Leg section – 90 -0 deg with gas spring Headrest -40 —+30 deg gas spring assisted Operational area around 500X1950 Height Capable to rise up to 1.2 meter. supplied with accessory like I.V rod, Lithotomy Legs etc.

### **Specifications :**

- Easy installation
- Accurate dimensions
- Optimum performance
- M.S. tubular frame work on 10 cms castors
- Height adjustment by crank mechanism with lever handle
- Five section S.S. tops
- All positions viz, "straight", Trendelenburg & Reverse Trendelenburg Positions, Kidney bridge, Lithotomy position
- Overall size: 72" x 20" x 32"
- Fine finished

---

labycare

www.labycare.com | Email: labycare@gmail.com