

Labycare-Crash Cart (SKSS834)



Hospital Crash Cart manufacturer by Labycare in India.

Description :

Labycare Hospital Crash Cart is a set of shelves on wheels used in hospitals for transportation and dispensing of emergency equipment at site of medical emergency for life support protocols to potentially save someone's life. The cart carries instruments for cardiopulmonary resuscitation and other medical supplies while also functioning as a support litter for the patient.

Specifications :

- Overall approx. size: 960 L x 500 W x 1545 H MM
- Frame work made up of S.S. tubes
- Two light weight polystyrene boxes each with three drawers
- Upper drawer with medicine container of different sizes
- Provision to hold oxygen cylinder made up of M.S.
- Massage board for cardiac patient are provided
- Six number colour hand out small bins at eye level
- Total three S.S. shelves with railings to carry medical equipment
- Provided with S.S. I.V. Rod, corner buffer
- The cart mounted on 125 mm dia castors 2 with brakes
- Also available in Epoxy Powder Coated Framework & S.S. Top

labycare

www.labycare.com | Email: labycare@gmail.com