

Labycare-Hospital Full Fowler Bed (SKSSFB809A)



Hospital Full Fowler Beds manufacturer by **Labycare** in India.

Description :

A hospital Full Fowler bed is designed to provide maximum positioning flexibility for patient comfort and clinical care. It allows adjustment of backrest, knee rest, and leg rest positions, supporting sitting and reclining postures. Commonly used in hospitals, it aids recovery, nursing care, and patient mobility.

Specifications :

Bed type: Full Fowler hospital bed

Operation: Manual / mechanical crank operated

Adjustable sections: Backrest, knee rest, leg rest

Frame: Mild steel tubular construction

Bed surface: Perforated CRCA sheet

Finish: Epoxy powder coated

Side rails: Optional / collapsible

Mattress size: Approx. 1980 x 900 mm

Load capacity: Approx. 150–180 kg

Application: General wards and nursing care units

labycare

www.labycare.com | Email: labycare@gmail.com

