

Labycare-Fume Hood (SKSS396A)



We have highly acknowledged organization engaged in presenting a remarkable range of Laboratory Fume Hood.

Description :

The window is designed into automatic balance type and is equipped with proper poise for taking its place easily. The bottom is equipped with stainless steel water tank with which liquid antiseptic, experiment leftover can be discharged from this tank through bathing so as to keep experiment environment safe and reliable. Thus, it is a kind of common lab equipment for discharging bad gas, cleaning and discharging in the course of experiment.

Laboratory furniture ventilation fume hood

Specifications :

Model NO.

SKSS396A

Customized

Customized

Condition

New

Place of Origin

INDIA (AMBALA)

After-Sales Service Provided

No Overseas Service Provided

Width

1.5m

Machine Type

Fume Hood

General Use

Commercial Furniture

Cerification

ISO 9001

Size

1500*760*1950mm

Trademark

LABYCARE

Transport Package

Wooden

Specification

1500 600 1800mm

Origin

INDIA

HS Code

94031000

labycare

www.labycare.com | Email: labycare@gmail.com