

Labycare-Slide warming table (SKSS350)



Slide Warming Table manufacturer by **Labycare** .

Description :

Slide Warming Table is very useful for the heating slides to a uniform constant temperature below the melting point of the wax. The surface temp is controlled by hydraulic type thermostat. Supplied complete with pilot lamp, cord, and plug but without thermometer.

Specifications :

- Works on AC 230 Volts 1 phase 50 Hz.
- Size: 150 X 600 X 150 mm
- Latest Model with complete body made of Stainless steel (304)
- Thermostatic Temperature controller
- Temperature Heating Indicator shows status of heating. Aluminum Top for best Uniform Heating. Gap of wall filled with glass fiber wool.
- **Other Specifications:**

Temperature Range	70° C
Length	12 inch
Width	10 inch
Height	8 inch

labycare

www.labycare.com | Email: labycare@gmail.com