

Labycare-Table top centrifuge (SKSS353)



Table top centrifuge manufactured by **Labycare** in India. Innovative designs enable safe, simple and comfortable operation every day.

Description :

Labycare Table top centrifuge is an Ideal for routine processing, the new ST 40 centrifuge series combines exceptional capacity, ergonomics and value with marketing-leading capabilities including unparalleled productivity and throughout as well as exceptional ergonomics. The Sorvall ST 40 series can process over 40% more 50 mL conical tubes, or 100% more blood tubes, with biocontainment sealing, than competitive models. Innovative designs enable safe, simple and comfortable operation every day.

<u>Specifications :</u>	Size/Dimension	390 x 470 x 280(H)mm
	Brand	Labycare
	Model Name/Number	SKSS353
	Capacity	300 ml, 50 ml x 6 Tubes(Swing
	Material	ABS
	Voltage	12V
	Phase	Single Phase
	Frequency	50 Hz
	Timer	99 min, Hold for long run
	Weight	20 Kg
	Max Rcf	2,790 x g(Swing Rotor) / 2,415
	Standard Rotor	Swing Rotor VR-45S
	Operating Temperature	-5 Degrees Celsius to +25 Deg
