

Labycare-VDRL shaker (SKSS358)



VDRL Shaker manufacturer by **Labycare** in India. Having a rubber covered platform for keeping slides, blood bottles, flasks, etc. by using spring holder.

Description :

Main indicator cum on / off switch and **Digital RPM meter** from 150 to 350 RPM with manual timer. Compact shaker for mixing blood samples and Wide speed range Timer options for accurate mixing time. It is also Digital RPM controller and digital timer available for exact results. Stable at high speeds Separate platform attachment available for mixing more than one samples.

Specifications :		
Type of Movement		Rotational
Speed Range		0 – 350 RPM
Operation Mode		Continuous operation
Timer Option		Yes
Dimensions		37 x 37 x 18 ½
Weight		11 kgs
Permissible Ambient Temperature		5 – 40 deg.
Permissible Relative Moisture		30%
Motor Type		PMDC motor
Maximum Load		6^ kg
Speed Display		None
Supply		220 V, 50 Hz.

