

Labycare-Laboratory Stirrer Tissue Homogenizer (SKSS360)



Laboratory Stirrer Tissue Homogenizer Manufacturer by **Labycare** in India.

Description :

1. Small in size, can save the shipping cost and takes a little space.
2. The small lab homogenizer is portable.
3. Though simple and small, the lab homogenizer can be used to mix various kinds of liquid, cream and paste material. It is also workable for high viscosity material.
4. The rotation speed is easy to be adjusted.
5. Height of the mixing rod/paddle is adjustable within a certain range.
6. Various homogenizer heads optional.

Specifications :

- **Application:** For Special all type of solvents used in Laboratories
- **Capacity:** Up 50 Litres
- **Type of Blade:** Special as per Requirement
- **Motor Hp:** Upto 2.0 Hp/ Available in Flameproof & Non flameproof version
- **Extra Accessories:** VFD, Control panel, Timer, Rpm Controller/Available in Flameproof & Non flameproof version
- **MOC:** Contact Parts- SS304, SS316, SS 316 L Alloy 20, Hastelloy C Non Contact part-SS304
- **Special Features:** Heavy Platform, UP & Down Motion with Gear box, RPM Controller, These type of Assemblies Can be Available in Both Flame Proof and Non Flameproof version

Documents : images/370382148ef2ff3e8aee6ab11eb980a8.jpg

labycare

www.labycare.com | Email: labycare@gmail.com