

Labycare-Environmental Chamber (SKSS367)



Environmental Chamber Manufacturer by **Labycare** in India. Being one of the trusted firms in the industry, we are highly engaged in providing a unique range of Environmental Chamber.

Description :

Labycare Environmental Chamber is also known as a climatic chamber. It is a system that enables users to manipulate the environmental conditions of an enclosed space to run controlled tests on a subject.

Specifications :

Technical Specification

- Outer body material of environmental chamber: Stainless Steel (304)
- Insulation Material: Asbestos
- Number of Rack: NA
- Humidity Chamber Capacity (in ltr): 300
- Thickness of insulated material (in mm): 35 Capacity of in built water tank (in ltr) 3
- Thickness of Inner Stainless Sheet (in mm): 0.9
- Thickness of outside pre-painted zinc coated sheet (in mm): 0.9
- Thickness of outside stainless: 1.0
- Four points temperature/ humidity adjustment: Yes
- Calibration-three freely selectable temperature value: Yes
- Fan-Control- depending on operation status automatic adaption of fan speed: Yes

labycare

www.labycare.com | Email: labycare@gmail.com

