

Labycare-Lypholizer (SKSS371)



Lyophilizer Manufactured by **Labycare** in India.

Description :

Labycare Lyophilizer also known as Industrial Lyophilizer. It is the process of removing the moisture from a biological product while maintaining the Integrity of the products biological & chemical structure and activities. WEIBER High Capacity Freeze Drier are conceptualized and designed to provide moisture removing solutions to all the microbiology, biochemistry, biotechnology and genetic engineering laboratories which require complete moisture removal and storage facilities of various biological agents/material. Beside the normal drying applications our freeze driers are fit to be used for various applications tests involving micro-organisms, plants, tissues etc. It is widely used in the medical, pharmaceutical, biological, chemical and food and other fields, FD series is suitable for laboratory use or low-volume production to meet the requirements of most laboratories conventional lyophilization.

<u>Specifications :</u> Power Consumption	230 V/50 Hz
Production Capacity	100 Kg
Brand	Labycare
Defrosting Method	Auto
Display	Touch
Drum	SS 304
Usage/Application	Laboratory
Voltage	230
Plate Load capacity	1 L
Vacuum Degree	10 Pa
Heating Power	950 W

