

Labycare- PLASMA FREEZER -80 (SKSSPF385)



Keyboard lock and password protected configuration page, Delayed start and safe stop interval between restart and being terminated Audible and visual alarm: High or low temperature alarm, Power failure alarm, Low battery alarm, Door open alarm, Filter blocking alarm, system failure alarm. With temperature printer, record 7days data. 72 hour battery backup for printer and controller display. Power supply: 220V /50Hz 1 phase, can be changed as 220V 60HZ or 110V 50/60HZ.

Description :

Keyboard lock and password protected configuration page, Delayed start and safe stop interval between restart and being terminated Audible and visual alarm: High or low temperature alarm, Power failure alarm, Low battery alarm, Door open alarm, Filter blocking alarm, system failure alarm. With temperature printer, record 7days data. 72 hour battery backup for printer and controller display. Power supply: 220V /50Hz 1 phase, can be changed as 220V 60HZ or 110V 50/60HZ.

Specifications :

- . Micro Controlled Temperature controller with Audio Visual
- . Alarm for high and Low Temperature
- . Heavy Duty Stainless Steel (304) Basket.
- . Display of Set and Process value.
- . Extra-strength non CFC refrigeration System.
- . The unit is mounted on castor wheel for easy mobility.

. CFC Free PUF (Poly Urethane Foam) insulation to ensure temperature stability and save power consumption.

labycare

www.labycare.com | Email: labycare@gmail.com