

Labycare-Cadavers Tank (SKSSCT380)



Cadavers Tank

Description :

- Strong Stainless Steel frame work construction. All structure polished finishing.
- Tray constructed with heavy duty stainless steel material for heavy drainage of the chemical.
- Provision of upside open cover.
- Tank Made from Thick Stainless Steel, with a high boarder.
- Provision of both side manual lifting handles.
- Provision of lock tray at the bottom.

Specifications :

Brand	SKSS
Type	Manual
Material	SS
Size	600 mm
Dia Height	900 mm
No of Wheel	4

labycare

www.labycare.com | Email: labycare@gmail.com