

Labycare-Dissection Table (SKSSDT383)



Dissection Table

Description :

Dissection Table

-

Specifications :

- Frame work constructed with stainless steel material.
- Top is constructed with stainless steel sheet, with high boarder.
- Provision for drainer hole with tube for waste liquid.
- Provision for Keeping SS bucket under the Top surface.
- Over all size: - 72"(L) x 24"(W) x 90(H) mm.
- Mounted on 125 mm swivel castors, two equipped with breaks.
- Supplied in knock down position, with assembly Kit.
-

labycare

www.labycare.com | Email: labycare@gmail.com