

Labycare-Plasma Expressor Mannual (SKSSPEM390)



SKSS is the largest manufacturer and exporter of Plasma Expressor in India.

Description :

SKSS is the largest manufacturer and exporter of Plasma Expressor in India.

Specifications :

- . Large plate supports bags up to 600 ml.
- Audible alarm.
- Indicator lamps.
- Operate 5 electronic expressors on one A/C power adaptor.

- Adjustable optical sensor detects red cells.
- Electronic clamp activates when optical sensor detects red cells.
- Enables operator to multi-task during plasma separation process.
- Provides uniform production output.

- Available in Stainless steel and powder coated mild steel cabinet withstand corrosion, wear & tear on continuous usage.

- The Compression plate is designed to exert uniform pressure on the blood bag.

labycare

www.labycare.com | Email: labycare@gmail.com