

## Labycare-Cryo Bath (SKSSCB391)



Cryo Bath

### **Description :**

Cryo Bath

### **Specifications :**

- Audio visual alarm if temperature deviates from preset temperature with silence function key with LED display.
- 0 Power fail indication with alarm.
- Prob fail indication with alarm.
- Display of Set value & process value & High - Low temp setting with visual alarms.
- Smooth acrylic tray for accommodating Plasma Bags.
- Rounded corners for better water circulation, improved uniformity and easier cleaning.

- Foamed in place. PUP insulation improves temperature stability and energy efficiency.
- 0316 stainless steel construction resists chemicals, simplifies cleaning.
- Long-life heating elements for quick response to controller.
- Attractive front panel bezel for easy cleaning.
- Cryo Bath size 4/6/8 Bags as per customer requirement.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)