

Labycare-Plasma Bath (SKSSPB392)



Plasma Bath

Description :

- Microprocessor based precise temperature control
- High performance 1.5 kw heater & compressor to achieve the set temp
- Over temperature safety limit, audio visual alarm
- LED monitor to display the working status of thawing bath
- Continuous water circulation for temp uniformity & rapid thawing
- Audio visual alarm after elapsed time

Specifications : Temperature	+ 37°C (accuracy $\pm 0.2^{\circ}\text{C}$)
External Dimensions W x D x H (mm)	760 x 530 x 640
Thawing Timer	Adjustable up to 2 Hrs
Internal Body Material	Stainless Steel AISI304 Grade
External Body Material	Powder Coated CRCA Steel
Noise Level	Less Than 65 db(A)
Recommended Line Voltage Corrector	3 KVA
Electrical	220-240 Volts, 50 Hz, SINGLE PHASE

labycare

www.labycare.com | Email: labycare@gmail.com