

Labycare-Dispensing Both (SKSS397)



Dispensing Both manufacturer by **Labycare** in India.

Description :

_____ is generally used by all pharmaceutical companies when raw material needs to be sampled and dispensed within a clean environment. They are designed to provide operator and products protected from outdoor or external contamination.

Specifications :

- Cleanliness - Class 100
- Particle Retention- 0.3 micron
- Velocity- 100fpm + 20%
- Illumination- 4 ft. 40w with Diffuser
- Noise Level- 60-65
- Decibels Standard- fed 209e
- Power Supply - 220v single phase
- Mangnehile Gorge
- UV light
- Hour Meter
- Brand- SKSS
- Model Number- SKSS397

labycare

www.labycare.com | Email: labycare@gmail.com