

Labycare-Semi Fowler Bed With S.S. Panel (SKSS809L)



Semi Fowler Bed With S.S. Panel manufacturer by **Labycare** in India.

Description :

Labycare Semi Fowler Beds with With S.S. Panel is utilized for backrest functions and knee rest. These fowler beds are furnished with collapsible side rails or simple lifting to guarantee the patient security. Semi Fowler beds are made up of Rectangular Tube/Round Tube. Sk Scientific & Surgical is outstanding Supplier giving different varieties of Hospital furniture like Semi Fowler beds everywhere throughout the world at discount aggressive costs from India. The different Positions can be taken by Smooth Crank System with fordable handles. To match up the changing and growing needs of customers, Sk Scientific & Surgical are engaged in offering a large range of Semi Fowler Bed with S.S. Panel. This Semi Fowler Bed with S.S. panel provided by us is used in hospitals and many more places. We offer product to the customers as per their needs and quality standards. Moreover, customers can avail customized product from us by giving their specifications on request. Sk Scientific & Surgical is considered as a renowned manufacturer and supplier of Hospital Bed.

Specifications :

- ABS moulded Head and Foot Panel with locking and corner buffer
- Framework made up of M.S. tubes
- S.S. Head & Foot Panels

- 2 section top made up of oval perforated CRCA M.S. sheet
- Legs fitted with PVC stumps
- Four location for I.V. Rod
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2000L x 910W x 600H mm
- Positions obtained by Smooth Crank System with foldable handle:

1. Backrest elevation

Optional

- 2 Section Mattress
- S.S. Collapsible Railing
- S.S. Telescopic I.V. Rod
- Advance linear actuator system (01 Pc Motor, 01 Pc control box, 01 Handset)
- M.S. or S.S. Head & Foot Panels
- Overall size 198L x 90W X 60H cms
- Finish: epoxy powder coated

labycare

www.labycare.com | Email: labycare@gmail.com