

Labycare-Pediatric bed (SKSS809H)



Pediatric Bed manufacturer by **Labycare** in India.

Description :

Labycare Pediatric Bed are especially designed for the Paediatric and Neonatology Department of hospitals, nursing homes and other healthcare units. Our Pediatric Beds ensure complete safety of the baby.

Specifications :

- Tubular Head and Foot Bows of unequal height with vertical tube support
- Frame work made up of M.S. tubes
- Top made up of oval perforated CRCA M.S. sheet
- Legs fitted with PVC stumps
- Four location for I.V. Rod
- Full Length drop side rails on both sides
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 1460L x 830W x 600H mm
- Bed frame approx size: 1370L x 760W mm

Optional

- Plain Mattress
- S.S. Telescopic I.V. Rod
- Full Length S.S. Side Rails

labycare

www.labycare.com | Email: labycare@gmail.com