

## Labycare-Labour bed with adjustable function (SKSS819)



**Labour Bed with adjustable function** manufacturer by **Labycare** in India.

### Description :

**Labour Bed with adjustable function** has been designed to provide comfort and safety to the mother during the crucial time of child birth. It is a highly flexible examination couch that is readily converted into a bed. It is designed for all stages of labour, delivery and recovery. Detachable polymer Moulded head and foot bows provide a soothing look. Provided with polymer Moulded safety side rails. Corner PU buffers absorb impact and reduce damage in handling. Provided with heavy duty saline stand with provision for fixing it at all four locations. Castors of 125/150 mm dia. With central locking mechanism for mobility with control. Back rest, Trendelenburg/reverse Trendelenburg, high – low, position maneuverable by separate actuator operated by patient control. The complete metal components are 7 tank pre treated and epoxy powder coated. Provided with high density three section mattress. Accessories: one collection tray, one set of lithotomy rod and crutches, one heavy duty in pole, one set of hand grip.

### Specifications :

- Overall Size: 1800L x 660W x 600-800H mm
- Tubular frame work mounted on 125 mm dia central locking castors
- Three section top made up of Bakelite Sheet with 3 section mattress
- Trendelenburg / Reverse Trendelenburg positions
- Adjustment through gas spring
- Height adjustment through hydraulic pump

- Back section adjustment through crank mechanism
  - Leg end section can slide under the main section
  - Provided with a pair of Upholstery Lithotomy positions
  - S.S. Telescopic IV Rod
  - Detachable S.S. Tray
- 

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)