

Labycare-Delivery Bed (SKSS818B)



Delivery Bed manufacturer by **Labycare** in India.

Description :

We are involved in offering best quality of **Labycare Delivery Bed**. The offered Delivery Bed is extensively used in hospitals. Further, the offered products are available in market in different sizes and specifications. Moreover, we ensure to deliver these products in various places.

Specifications :

- Height adjustable by a pair of knee crutches
- Approx. dimensions: 72" L x 30" W x 30" H
- Tubular frame work made in two sections
- Body section with legs mounted on PVC stumps
- Leg section on swivel castors with brakes
- Two section cushioned top
- Finish: Pre-treated & Epoxy Powder Coated

labycare

www.labycare.com | Email: labycare@gmail.com