

Labycare-Bed Side Screen (SKSS851A)



Bed Side Screen manufacturer by **Labycare** in India.

Description :

Labycare Bed Side Screen is used in operation theaters as well as in dispensaries, medical clinics and medical shops. This helps the doctors and chemists to maintain privacy for their patients.

Specifications :

- Overall size: 1680H x 2450L mm
- M.S. tubular frame work made up of 19mm x 18G in four section
- Legs made of 25mm x 18G mounted on 5 cms castors with hooks and spring along-with curtains
- Each span is 610 mm wide
- Finish: Pre-treated and Epoxy Powder Coated

labycare

www.labycare.com | Email: labycare@gmail.com