

Labycare-Carm OT Table (SKSS854)



Carm OT Table manufacturer by **Labycare** in India.

Description :

Labycare Carm OT Table is five section table top gives ease for various surgical positions. Its _____

Specifications :

Technical Specification

- Overall Length of the table :- 1950mm
- Overall width of the table :- 530mm
- Minimum height of the table :- 780mm
- Maximum height of the table :- 1030mm
- TB/RTB :- 40 degree
- Flex/Reflex :- 80/220 degree
- Side Tilt :- 22 degree
- Kidney section :- 125 mm manually
- Leg Plate :- 90 degree Down
- Hi-Low Lift :- Electric

Mechanical Movement

- Flex-Reflex
- Lateral tilt
- Trendlenburg /Reverse Trendlenburg

- Leg rest
- Head rest
- Kidney bridge
- Floor lock

labycare

www.labycare.com | Email: labycare@gmail.com