

## Labycare-Hydraulic OT Table (SKSS855 )



**Hydraulic OT Table** manufacturer by **Labycare** in India.

### Description :

**Labycare Hydraulic OT Table** are designed to provide accurate positioning with an easy mechanical operation. These tables are ideal for use in all the major surgical process that requires precision.

<b>Specifications :</b> Material	Stainless Steel/ Mild Steel/ Cus
Brand	Labycare
Usage/Application	Hospital
Minimum Height	750 mm
Load Capacity	200 Kg
Battery Backup	Upto 3 Hour
Is It Adjustable	Adjustable
Electric Control	Yes
Width of table top	500 mm
Up-down movement	772-950 mm
Reverse Trendelenburg	25 Degree
Lateral tilt	20 Degree (Both Sides)
Leg section	90 Degree (Down Manual)

---

labycare

www.labycare.com | Email: labycare@gmail.com