

Labycare-Hydraulic Operating Table (SKSS964)



Hydraulic Operating Table manufacturer by **Labycare** in India. Quick and stable movements materials, durable and reliable Most convenient control switch and easy transfer.

Description :

The Operating Table is a crucial component in surgical settings, designed to provide a stable and adjustable surface for the patient during surgical procedures. This Universal OP Table offers a combination of mechanical and hydraulic features for versatile positioning and ease of use.

Specifications :

Technical Specifications

- **Type:** Hydraulic Operation Table
- **Material:** Stainless Steel (SS 304)
- **Table Top:** X-ray translucent, multi-section
- **Height Adjustment:** 650 – 1040 mm
- **Dimensions:** Approx. 2000 (L) × 500 (W) mm
- **Operation:** Foot pedal (height), manual controls
- **Lateral Tilt:** ±20°
- **Trendelenburg:** Up to 45°
- **Reverse Trendelenburg:** Up to 20°
- **Back Section Adjustment:** +75° / -45°
- **Head Section Adjustment:** +50° / -90°
- **Leg Section Adjustment:** +22° / -90°
- **Longitudinal Shift:** 0 – 300 mm (for C-arm compatibility)
- **Safe Working Load:** 200 – 250 kg
- **Base Type:** Heavy-duty mobile base with locking system
- **Mattress:** High-density foam, waterproof & anti-static cover
- **Finish:** Matte / satin finish, corrosion resistant
- **Castors:** Anti-static, smooth rolling with brakes
- **Imaging Compatibility:** Suitable for C-arm / X-ray procedures

Optional Accessories / Configurations

- Anesthesia screen frame
- Shoulder supports
- Arm boards (adjustable)
- Lateral supports
- Kidney bridge attachment
- Lithotomy leg holders
- Orthopedic traction attachments
- Headrest (neurosurgery / ENT type)
- Gel pads / pressure-relief mattress
- Stainless steel instrument tray
- Central locking system upgrade
- Battery backup (for hybrid models)

labycare

www.labycare.com | Email: labycare@gmail.com