

Labycare-Shadow Less Light (SKSS858)



Medical multifunctional LED Operating Shadowless Lamp Hospital. a long time for work, and does not emit excess heat, because of overheating will cause the operation of discomfort, so as to enhance the probability of successful operation of the doctor.

Description :

Medical multifunctional LED Operating Shadowless Lamp Hospital. a long time for work, and does not emit excess heat, because of overheating will cause the operation of discomfort, so as to enhance the probability of successful operation of the doctor.

Specifications :

1. The lamp cap adopts the fully enclosed and non-angular streamline design and unique spraying technology to ensure the surface of the lamp cap is corrosion-resistant, anti-collision, easy to clean, no screw exposed, and good laminar flow effect.
2. the unique polygonal mirror design, using the lens principle of first reflection and refraction, so that the light source of each reflective system is just rightly collected into the operation area, even if the entire lamp head is completely blocked by the doctor's head, it can reach 100% without shadow. effect. Avoid the glare and light decay caused by direct light, and avoid the damage and visual fatigue caused by the direct light source to the doctor's eyes.
3. high-performance LED light engine, modular design, easy to replace.
4. with minimally invasive surgery mode one-button switching, providing illumination during endoscopic and laparoscopic surgery.

5. the lamp head operation is light and flexible, accurate positioning, no drift.

labycare

www.labycare.com | Email: labycare@gmail.com