

## Labycare-Orthopedic operating table (SKSS240)



**Orthopedic Operating Table** manufacturer by **Labycare** in India.

### Description :

**Labycare Orthopedic Operating Table** are designed for orthopedics operation and surgical procedures according to clinical requirements of orthopedics surgeries. A patient lies on the operating table during a surgical procedure. The purpose of a surgical table is to keep the patient in place while the surgical team operates, and may move various parts of the body using surgical table accessories for easier access to the surgical site.

### Specifications :

- Type of OT table Actuator : Electrical
- Table top length, Range : 1900 mm-2100 mm
- Table top Height, Range : 650 mm - 1150 mm
  
- Provision of X ray cassette channel : Yes
- Table top longitudinal sliding, Range : 5" to 12"
- Kidney Bridge required : Yes
- Split leg section , detachable-1 pair : Yes
- Battery backup of operating table in hrs : 3 or more than 3 hrs
- Battery backup provided : Yes
- Foot pump : No
- Extension bar-1pair : No
- Horizontal Attachments-1no : Yes
- Remote-1 no : Yes
- Orthopedic Attachments : Floor based
- Foot plate mount : 2
- Supports bars : 2
  
- Brand : SKSS

- Model No : SKSS240

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)