

Labycare-Deep Freezer For Universal Immunization Program (SKSSDFFUIP101)



Deep Freezer For Universal Immunization Program manufacturer by **Labycare** in India.

Description :

Labycare Deep Freezer For Universal Immunization Program is used to maintain temperature. This is used for storing of OPV (district level and above only) and also for freezing ice packs.

Specifications :

- **Warranty of Deep Freezer:** The Equipment should have minimum warranty including comprehensive maintenance of sixty months after installation or sixty months after the supply which ever is late
- **Lower Limit Temperature Range:** -30° Celcius
- **Upper Limit Temperature Range:** +50° Celcius
- **Temperature Accuracy:** +/- degree Celcius
- **Warranty of Alcohol Stem Thermometer:** The warranty shall remain valid for two years from the date of final acceptance
- **Automatic Voltage Stabilizer features:** Over load protection at input, automatic cut-off time between 2-4 seconds of Under voltage/Over voltage cut off at output ,Maximum power loss due to leakage of current should be less than 2% ,Restart delay of 6-9 minutes after the cut-off along with quick start button ,The output voltage should be stable and capable of operating continuously with respect to any change in ambient temperature and regulation should be automatic ,Transformers copper wound on CRGO laminations, as per ISO specification for high efficiency and low loss ,Voltmeter (digital/analog)with selector switch to indicate input and output voltage with off-position ,The output side shall be provided With C series MCB with braking current of 10KA duly ISI marked ,Provision for two output sockets of 15Amps ISI marked ,Details of transformer should be furnished,such as rating, weight of winding and core, which should be verified during inspection ,A Metal Oxide Varistor (MOV) for surge suppression should be provided and connected directly on the power line

labycare

www.labycare.com | Email: labycare@gmail.com