

## Labycare-Emergency Recovery Trolley (SKSS835)



**Emergency Recovery Trolley** manufacturer by **Labycare** in India.

### Description :

**Labycare Emergency Recovery Trolley** is designed in compliance with international standards. The product is designed and manufactured to give provide safety to patient and care giver. ERT has SS make drop down side rails and patient safety strap over the stretcher to lock the patient movement.

### Specifications :

- Two section top height adjusted by foot operated imported hydraulic pump from 660 mm to 910 mm
- Gas spring assisted Trendelenburg/Reverse Trendelenburg positions
- Complete with corner buffers, synthetic rubber covered handles, sliding X-ray cassette, accessories tray, oxygen cylinder holder
- S.S. Telescopic IV rod and swing away S.S side rails
- Finish: Pre-treated and Epoxy Powder Coated
- Four swivel castors, 5" dia, two with total lock
- X-ray translucent removable stretcher, backrest raised on ratchet
- Overall approx size: 210 L x 65 W x 65-95 H cms

---

labycare

www.labycare.com | Email: labycare@gmail.com