

Labycare-Mortuary Chamber (SKSS376)



Mortuary Chamber Manufacturer by **Labycare** in India.

Description :

Labycare Mortuary Chambers are low temperature refrigerated cabinets that are used to keep dead bodies for short or long time. These units are widely used at hospitals, railways, airports, disaster camps and army etc. These cabinets are also called by other names such as morgue freezer, dead body storage refrigerator and mortuary fridge. Usually the temperature of mortuary cabinet remains between 2°C to 5°C but in some cases freezing temperature up to 0°C, -10°C and -15°C can also be designed.

Specifications :

Temperature range	-15°C to 10°C
Temperature controller	<ul style="list-style-type: none"> • Microprocessor controller • High resolution LED display • Elegant flat keypad • IP65 protection • Made in india
Display precision	0.1°C
Temperature sensor	Platinum resistance temperature sensor
Overall construction	Stainless steel 304
Insulation	PUF

Trolley

- Made of stainless steel
- Run on rollers fitted with ball bearings

•

	<ul style="list-style-type: none"> • Fixed on SS hinges • Spring door latch
Refrigeration system	<p>Individual refrigeration system for each chamber</p> <ul style="list-style-type: none"> • Karloshkar make compressor • Fluoride-free refrigeration system • R134a refrigerant • EBM Motor fan
Air circulation	Axial fan
Illumination	LED light
Alarms	Audio / visual alarms
Power supply	220 Volts 50Hz
Optional	<ul style="list-style-type: none"> • Battery backup for temperature & alarm • Temperature chart recorder • Voltage stabilizer

STANDARD SIZES / DIMENSIONS

MODEL	BODIES	COMPRESSOR UNITS	EXTERNAL DIMENSIONS
SKSS37 6a	ONE BODY	1	As per customer requirement
SKSS37 6b	TWO BODY	1	As per customer requirement
SKSS37 6c	THREE BODY	1	As per customer requirement
SKSS37 6d	FOUR BODY	1	As per customer requirement
SKSS37 6e	SIX BODY	1	As per customer requirement
SKSS37 6f	NINE BODY	2	As per customer requirement

labycare

www.labycare.com | Email: labycare@gmail.com