

Labycare-ROTARY MICROTOME (SKSSRM951)



ROTARY MICROTOME /SKSSRM951

Description :

ROTARY MICROTOME

Specifications :

Section Thickness Setting	From 0.5 μ m-2 μ m in 0.5 μ m increments From 2 μ m-20 μ m in 1 μ m increments From 20 μ m-60 μ m in 5 μ m increments
Trimmming Thickness	10 μ m-6 μ m
Vertical Stroke	60 m
Specimen Orientation	Horizontal 8°, Vertical 8°
Object Feed	28 mm
Max. Sectioning Area	60 mm x 50 mm
Dimensions	
Widht(Incl.Handwheel)	460 mm
Depth(Incl.Waste tray)	500 mm
Height(Total)	300 mm
Working Height(Knife Edge)	105 mm
Weight(without acc.)	23 kg. (Approx)

labycare

www.labycare.com | Email: labycare@gmail.com