

Labycare-CRYOSTAT MICROTOME (SKSS952)



CRYOSTAT MICROTOME

Description :

CRYOSTAT MICROTOME

<u>Specifications :</u>	Usage/Application	Hospital,Clinic
	Model	SKSS952
	Cooling Temperature	-40 degree C
	Temperature Range	-10C to -40C
	Speed	0.7mm/s,0.35 mm/s
	Weight	130 kg
	Frequency	50 hz
	Voltage	220v
	Automation Grade	Semi-Automatic
	Operating Voltage	220v,50Hz
	Warranty	12 months
	Dimensions	660x640x1130 mm
	Vertical Stroke	60mm
	Maximum Specimen Size	45x55 mm
	Specimen Orientation X And Y Axis	12 degree
	Nominal Voltage	220 v
	Specimen Orientation Z Axis	360 degree
	Thickness	1-100 um
	Trimming Thickness Range	10-400 um adjustable
	Is It Portable	Non Portable
	Max Specimen Size	45x55 mm
	Brand	Radical
	Glass Viewing Window	yes

labycare

www.labycare.com | Email: labycare@gmail.com