

Labycare-Dissection Microscope (SKSSDM956)



Dissection Microscope manufacturer by **Labycare** in India.

Description :

Labycare Dissecting Microscope is used to view three-dimensional objects and larger specimens, with a maximum magnification of 100x. This type of microscope might be used to study external features on an object or to examine structures not easily mounted onto flat slides.

Specifications :

• Eye Piece	Pair WF 10x20 mm eye pieces.
• Head	Trinocular heads with 45° inclined Both eye pieces with 5 mm adjustments. Inter pupillary distance adjustable from 55 mm to 15 mm.
• Objective part	Zoom Ratio 6.4 : 1. 7x 10 45x magnification Field of view from 33 mm to 10.49 mm Working distance 100 mm.
• Stand	Ergonomically designed pillar Binocular stand with LED illumination clamps. Alloy metal cast, hardened off-white.
• Focusing Adjustment	Coarse adjustment, moving on Guide ways
• Illumination	Transmitted and incident 3W LED with built-in 100-240V power supply Both illuminators can be used separately and the light intensities can be adjusted separately.

labycare

www.labycare.com | Email: labycare@gmail.com